

## **FLEXBRITE II**

### **Bright White LED Illumination Adjustable Sharp Spot Light**

(Spot diameter: 130-300mm )  
40" Flex Arm


Dual reinforcement of flexible arm to avoid drifting  
Bright illumination of the examination area

**ETL approved for Canada & United States**  
Available with stainless Steel Wall mounting plate  
and weighted mobile base

MEL145W  
MEL145S

With Wall mount  
With Mobile base

Volume Discount available



Spot adjusting ring

On / Off Knob  
&  
Brightness adjustment



Weighted base

Lamp Power: 5W  
C.T.: 5500 ± 500 K  
Lamp Life: 50000 h  
Illumination: 1600 Lux  
Power supply: 110-130/220-240 V, 50/60 Hz