

# MFM-CMS

## Central Monitoring System

# EDAN

### Specification

<b>Max. Amount of Networked Monitors</b>		<b>Compatible Monitors</b>	
Up to 128 bedside monitors: Support 32 bedside monitors in one display Support 64 bedside monitors in two displays		elite V8, iM80, iM50, iM9, iM9A, iM8, iM8A, iM8B M80, M50, M9, M9A, M8, M8A, M8B M3, M3B	
<b>Alarm</b>		<b>Print</b>	
Alarm Categories:	Physiological/Technical	Patient Information	Waveform Review
3 Levels of Alarm:	High/Medium/Low Level	Alarm Wave	NIBP Review
Indicating Method:	Audible/Visual	Alarm Table	Drug Calculation, Titration Table
		Trend Graph	PDF File
		Trend Table	

### Recommended System Configurations

Component	Requirement	
PC Workstation	CPU: Memory: Hard Disk: USB Ports: Ethernet Card	Intel® Core™ 2 Duo or AMD Turion™ II M520 2 GB or above 320 GB or above 2 at least
Display	Resolution: LCD Display	4:3 / 16:9 (i.e. 1280×1024 / 1440×900)
Ethernet Router	Ethernet Ports: Cascade Supported	8 or above
Others	Keyboard Mouse Speaker Printer:	Support A4-size printing paper

### Ordering Information

<b>MFM-CMS</b>	LAN Solution base Central monitoring station Including standard accessories: MFM-CMS software CD, user manual Software installation manual USB sentinel, Quick reference card
<b>Wireless Kit</b> 1117047336 0201101206 0201101207	D-LINK wireless (DWL-2100AP) ASUS wireless AP(WL-330g EAP) for M50 ASUS wireless AP(WL-330g EAP) for M80



# MFM-CMS

## Central Monitoring System



**EDAN**  
Edan Instruments, Inc.



**Solic Medical Equipment**

Ph. 416-740-8885-Fax, 416-740-1883  
For Inquiries: sales@solicmedical.com



Care for Health