

22401

Basic Examination Table

Features:

- *For Basic Examinations and Treatments*
- *Steel frame*
- *Two drawers*
- *One shelf*
- *350 lbs capacity*
- *Strong 4-leg base in white colour*
- *padded 2-inch mattress*
- *Mattress is available in black colour only*



Specifications:

- **Height: 30"**
- **Width: 24.5"**
- **Length: 75"**

Shipping Dimensions:

- **Height: 16"**
- **Width: 26"**
- **Length: 78"**
- **Weight: 145 kg.**

***Table is shipped unassembled**