

# 32501

# Family Practice Exam Table

## Features:

- *Seamless premium vinyl upholstery*
- *Easy to clean - compliant to infection control*
- *Fully welded, powder-coated, all-steel frame*
- *Pneumatic backrest Leg extension*
- *Pull-out stirrups that folds/slides away when not in use*
- *Table paper holder accommodates 18" or 21" paper widths*
- *Non-Slip retractable footstep*
- *Adjustable feet for leveling*
- *Weight capacity: 500 lbs.*

## Specifications:

- **Length: 58.5"**
- **Width: 27.5"**
- **Height: 34"**
- **Length w/ foot extended: 75"**
- **Base width: 22.8"**

### **Step Stool Dimensions:**

- **Length: 16.5"**
- **Width: 18"**
- **Height: 9.5"**

### **Drawer Dimensions:**

- **Length: 20"**
- **Width: 16.8"**
- **Height: 4"**

