

# 8031090/YXH1A6S

## CPR Board and Spinal Board

### Features:

#### 8031090-CPR Board:

- *Compression plate composed of high-density polyethylene plastic*
- *One concave end to support patient's head in comfortable position*
- *Easily portable*

### Dimensions:

- **Length: 22" x Width: 16" x Thickness: 2"**
- **Net Weight: 3 k.g.**
- **Gross Weight: 4.4 k.g.**

### Shipping Dimensions:

- **Length: 24.2" x Width: 17.71" x Height: 6"**

#### YXH1A6S- Spinal Board

- *Compatible with most head mobilization devices and strap mechanisms*
- *Includes black straps*
- *Composed of strongly-engineered plastic*
- *Easily stored*
- *Load-bearing capacity: Approximately 150 k.g.*

### Dimensions:

- **Length: 72" x Width: 21.18" x Thickness: 2"**
- **Net Weight: 7.5 k.g.**
- **Gross Weight: 8.8 k.g.**

### Shipping Dimensions

- **Length: 72.5" x Width: 18.5" x Height: 2.7"**



8031090

YXH1A6S