

AG4H

Emergency Transport Stretcher

Features:

- *Made of high-strength aluminum alloy*
- *Single-person operative*
- *Maximum Angle of the Back: 65°*
- *Foldable legs which are controlled by handles located on each end of the stretcher*
- *Weight capacity of 159kg.*

Specifications:

High Position:

- **Length: 74.8"**
- **Width: 21.5"**
- **Height: 36.2"**

Low Position:

- **Length: 74.8"**
- **Width: 21.5"**
- **Height: 9.8"**



Shipping Specifications:

- **Net Weight: 40kg.**
- **Gross Weight: 46kg.**
- **Carton Size: (1pc/carton): 195 x 62 x 25.5cm**