

EZ108

Foldable Stretcher

Features:

- *Made of high-strength aluminum alloy*
- *High strength Aluminum Alloy and Oxford Leather material*
- *Light weight*
- *Easily portable*
- *Flexible handles, and carrying bag included*
- *Weight capacity: 350 lbs.*

Specifications:

Unfolded:

- **Length: 76"**
- **Width: 21"**
- **Height: 4.2"**

Folded:

- **Length: 40"**
- **Width: 7"**
- **Height: 6"**
- **Net weight: 10.5 lbs.**

