

F-45-1 EMS Crash Cart

Features:

- 5 locking drawers
- Central lock system
- Includes CPR board and built-in telescopic IV pole
- One sharp cannister brackets
- 4 noiseless locking castors with breaks
- Slide-out stainless steel shelf, plastic top surface, and stainless steel top rail
- One oxygen tank holder

Specifications:

Cart Dimensions:

- Width: 32" x Length: 20.5" x Height: 39"

First & Second Drawers:

- Width: 20" x Length: 14.25" x Height: 2.5"

Third & Fourth Drawers:

- Width: 20" x Length: 14.25" x Height: 4"

Fifth Drawer (located at bottom of cart)

- Width: 20" x Length: 14.25" x Height: 8.5"

Shipping Dimensions:

- Length: 36.5" Width: 24" x Height: 36.5"
- Weight: 80 kg.

