

F-46 EMS Crash Cart

Features:

- 5 locking drawers
- Stainless steel top rail and top surface
- Plastic tray under the top surface
- Defibrillator holding tray (stainless steel with straps) and built-in I.V. pole
- Paper/chart holder
- 2 garbage cans
- 4 noiseless locking castors with breaks

Specifications:

Cart Dimensions:

- Width: 32" x Length: 20.5" x Height: 39"
- Defibrillator SS tray installed: 44.5"

First & Second Drawers:

- Width: 20" x Length: 14.25" x Height: 2.5"

Third & Fourth Drawers:

- Width: 20" x Length: 14.25" x Height: 4"

Fifth Drawer (located at bottom of cart)

- Width: 20" x Length: 14.25" x Height: 10.5"

Shipping Dimensions:

- Length: 37.5" x Width: 24.5" x Height: 40"
- Weight: 95 kg.

