

HX400C/HX200A/HX300

Height Measuring Equipment

Features:

HX400C- Tape Retractable Stadiometer:


- *Saves space with roll-up mechanism*
- *Easy installation*
- *Durable metal tape with length of 200cm.*

HX200A-Telescopic Stadiometer:

- Mechanical telescopic measuring rod with large measuring range
- Designed specifically for wall mounting
- The telescopic measuring rod ensures that the result can be easily read at eye level
- Measuring slide can be folded down for safety.
- 60 - 200cm

HX300- Infanometer:

- Most effective method of measuring a baby's growth
- Total measuring length: 8" - 39"



HX300

HX400C

HX200A