

J101A

Crash Cart, 5 Drawers

Turquoise/Gray

Features:

- *High-impact ABS material frame with four aluminum alloy columns*
- *Stainless steel guard rails*
- *Handles on both sides for convenient handling*
- *Two small drawers, two medium drawers, and one large drawer, all equipped with removable dividers for easy cleaning*
- *Dust bin (1 pc)*
- *ABS Side basket (1pc)*
- *Sharp container (1pc)*
- *Adjustable I.V. pole, 4 hooks (1 pc)*
- *Oxygen tank holder (1 pc)*
- *Sliding side shelf (1 pc)*
- *Transparent file box (1 pc)*
- *Two keys for centralized locking (1 set)*
- *Casters (5 inches), two with brakes*



Large Drawer:
52 x 36 x 12 cm
Panel Height: 21 cm

Middle Drawer:
52 x 36 x 12 cm
Panel Height: 14 cm

Small Drawer:
52 x 36 x 5 cm
Panel Height: 6 cm

Specifications:

- *33.4" x 20.4" x 38"*
- *Wooden Box packing: 36" x 23.4" x 40.7"*
- *Net Weight: 50 kg*
- *Gross Weight: 74 kg*