

LS-MIF

Fiber Optic Laryngoscope

Features:

- *Macintosh Fiber-optic Laryngoscope set*
- *High quality stainless steel construction*
- *Bright White light up to 10000 LUX*
- *Integrated fiber optic bundle*
- *Fully compliant to infection control guidelines*
- *Easy to clean and sterilize- No crevices or inaccessible areas*

Components:

- *1 x Handle- takes two C size batteries (not included)*
- *1 blade x 12232 Mac # 2 116mm*
- *1 blade x 12233 Mac # 3 135mm*
- *1 blade x 12234 Mac # 4 155mm*
- *1 blade x 12235 Mac # 5 178mm*
- *1 x User manual*



Health Canada approved
Distributed exclusively by Solic Medical in Canada