

MEL145

Flexbrite II Spotlight

Features:

- 40" Flex Arm
- Dual reinforcement of flexible arm to avoid drifting
- Bright illumination of the examination area
- Available with stainless steel wall mounting plate and weighted mobile base

Specifications:

- Spot diameter: 130-300mm
- Lamp Power: 5W
- C.T.: 5500 ± 500 K
- Lamp Life: 50000 h
- Illumination: 1600 Lux
- Power supply: 110-130/220-240 V, 50/60 Hz

Options:

- MEL145S: Includes mobile base
- MEL145W: Includes wall mount with clamp

