

NJ801

Crash Cart Deluxe, 5 Drawers Multicolour

Features:

- High-impact ABS material frame with aluminum alloy plastic plate, cold rolled steel frame, durable outer double-layer frame
- Handles on both sides for easy maneuvering
- Two small drawers, two medium drawers, and one large drawer, all equipped with removable dividers for easy cleaning
- Drawer with transparent medicine card
- Defibrillator shelf (1pc)
- Dust bin (2pcs)
- ABS Side basket
- Sharp container
- Adjustable I.V. pole, 4 hooks on top
- CPR Board at back
- Oxygen tank holder on side
- Sliding side shelf
- Two keys for centralized locking
- 5" Casters, two in front with brakes



Large Drawer:

52 x 36 x 12 cm Panel Height: 21 cm

Middle Drawer:

52 x 36 x 12 cm Panel Height: 14 cm

Small Drawer:

52 x 36 x 5 cm Panel Height: 6 cm

Specifications:

- 26.3" x 20.8" x 40"
- Wooden Box packing: 30" x 26" x 40.7"
- Net Weight: 56 KG
- Gross Weight: 78 KG