

AYR1009/AYR1009E

Three-Section Power Exam Table

Features:

- 17.5" profile in low position allows multiple-level transfers
- 80° repositioning attainable on both sides with pneumatic backrest
- Chair-to-Table design takes patient from consultation to examination level
- Heavy duty, all welded powder coated steel frame with 2" firm density foam tops
- Hands-free low voltage foot control
- Paper roll holder
- 500 lb weight capacity
- Built-in casters for easy movement

Specifications:

- Height: Max: 32"; Min: 17.5"
- Width: 27"
- Length: 78"
- Head Section: 39.5"
- Seat Section: 24.5"
- Foot Section: 13"

Options:

- **AYR1009:** Three-section table available without stirrups
- **AYR1009E:** Three-section table available with stirrups



AYR1009