

# Pocket spot checks at all times

## iM3s

Vital Signs Monitor



# Small outside, smart inside

Being the smallest vital signs monitor among its kind, iM3s presents a fashionable handheld design, delivering extraordinary portability throughout spot checks. Though small in size, it integrates all necessary vital measurements together, including SpO<sub>2</sub>, NIBP and optional TEMP.

## A Dedicated Spot Monitor

With its pleasing size and large capacitive touch screen, iM3s could be operated by one hand, bringing comfortable and familiar user experience to daily applications. For each test, all measurements and custom parameters could be displayed on the same screen, while the test results being saved automatically in its built-in non-volatile memory or transfer to hospital information system.



Lightweight Design



Waterproof & Dustproof



Shock Resistance



Minimum Run Time



Network Connection



Wireless Connection



Internal Storage



Barcode Scanner



# Tabletop Expansion

Plugged onto an optional tabletop dock, the iM3s easily fits into regular tabletop applications, enabling more possibilities, including AC charging, optional thermal recording, optional temporal temperature measurements, etc.



## Powerful functions with flexible configuration

### Main Unit

- Standard Configuration: SpO<sub>2</sub>, NIBP, Capacitive touch screen, Li-ion battery
- Optional Configuration: Built-in barcode scanner, 2.4G/5G Hz Wi-Fi, e-Link

### CS-04 Charger Stand

- AC power supply

### CS-05 Charger Stand

- Standard Configuration: Built-in Li-ion battery, AC power supply, USB port, RJ45 port
- Optional Configuration: Thermal recorder, Exergen TEMP (TAT-5000S)



## iM3s

Vital Signs Monitor

### About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.

 **Solic Medical**  
Ph: 416-740-8885  
info@solicmedical.com



  
**EDAN**

A world of potential