

# SD3

Ultrasonic Pocket Doppler



**Solic Medical**  
Passion to lead



A world of potential

# SD3

## Ultrasonic Pocket Doppler

With excellent sensitivity, complete interchangeability and high durability, SD3 series Ultrasonic Pocket Dopplers are ideal for routine fetal heart rate detection by clinicians.



### Configurations

Model	SD3 VASCULAR	SD3 LITE	SD3	SD3 PLUS	SD3 PRO
FHR Display	×	×	√	√	√
Charger Stand	×	×	×	○	○
Record & Play	×	×	×	×	√
AA Batteries	○	○	○	×	×
Li-ion Battery	×	×	×	√	√

√ Standard    × Not available    ○ Optional

## User-Friendly Design

- Compact design with easy one-hand operation
- Bright OLED screen with real-time FHR display
- Hi-Fi speaker

## Interchangeable Waterproof Probe


- Plug-and-play
- 2 MHz probe with deep penetration for later pregnancy and larger patients
- 3 MHz probe detects fetal heart beat as early as 9 weeks
- 4/5/8 MHz probe for blood flow detection
- Automatic power-off while putting the probe back to the rack (SD3 PLUS & SD3 PRO)

## Work Mode

- Real-time work mode for spot check
- Record & play with audio playback and numeric display

## Power Source

- Standard AA batteries for easy replacements
- Built-in Li-ion rechargeable battery ensures long working time
- Unique charger stand and wall mount design



SD3 VASCULAR



SD3 LITE



SD3



SD3 PLUS



SD3 PRO



## SD3

Ultrasonic Pocket Doppler

### About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**Solic Medical**

Ph: 416-740-8885  
info@solicmedical.com



A world of potential