

# A Slim Monitoring Experience

## iM3

Vital Signs Monitor



**Solic Medical**

Passion to lead



**EDAN**

A world of potential

# Smarter, Larger, yet Slimmer

Combining a stunning ultra-slim design, a large capacitive touch screen, together with a rich set of leading functions, iM3 takes the lead and presents itself as the beginning of a new generation of vital signs monitors. Bringing in full options of SpO<sub>2</sub>, NIBP and Quick Temp/Infrared Ear Temp and differentiated work modes, iM3 serves for both long time monitoring and spot check applications.



Touch Screen



- **Standard Parameters**

EDAN SpO<sub>2</sub>, EDAN NIBP, PR

- **Optional Parameters**

Nellcor™ SpO<sub>2</sub>, SunTech® NIBP, EDAN Quick TEMP, Infrared Ear TEMP, Covidien Quick TEMP

## Redefinition of Adaptive Work Modes

- Round Mode

Special work mode designed for daily ward rounding application. Worklist could be easily input through USB disk or downloaded from EMR system. Moreover, the color-coded patient list comes with square windows interface design, you can tell patients' status and information in a quick glance.



Patient management



Patient status with color code

- Monitoring Mode

For continuous patient monitoring with real-time data, real-time alarms, and history trend.

- Spot Mode

For spot-check applications, such as out-patient, emergency, etc. Real-time vitals are acquired and displayed, most of the alarms will be muted.

## TEMP Module Options



Infrared Ear TEMP



EDAN Quick TEMP



Covidien Quick TEMP



Rolling Stand

## Non-volatile built-in Memory

Spot Mode

16,000,000

Measurements

Round Mode

800,000 sets

Round Records

Monitoring Mode (for each patient file)

200

Alarm

1200

NIBP Measurements

240 h

Trend



## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**Solic Medical**

Ph: 416-740-8885  
info@solicmedical.com



**A world of potential**