

iM70 Patient Monitor



Solic Medical
Passion to lead



A world of potential

iM70

Patient Monitor



12.1"

Maximum 13 Waveforms



No-fan Design



Night Mode



Pitch Tone



Thermal Recorder



Shortcut Menu



Ethernet Printer



Standard Parameters

- 3/5-lead ECG, HR, RESP, EDAN SpO₂, EDAN NIBP, PR, 2-TEMP

Optional Parameters


- 12-lead ECG, 2-IBP(with ICP support), C.O., EDAN G2 CO₂ (Sidestream), Respiration CO₂ (Mainstream and Sidestream), MASIMO Anesthesia Gas/O₂, Nellcor OxiMax™ SpO₂, SunTech® NIBP



Reliable Algorithms

- iSEAP™ ECG algorithm optimized for arrhythmia detection, pacemaker detection, and HR measurement
- SEMIP® 12-lead diagnosis algorithm with 208 kinds of analysis results
- iMAT™ SpO₂ algorithm with outstanding motion resistance and low perfusion resistance performance
- iCUFS™ NIBP algorithm optimized for cardiac patients, hypertensive patients, and neonatal patients

Multiple Display Modes



Standard Display

Large Font

Vital Display

Trend View

OxyCRG

Bed View

End tidal Carbon Dioxide for Intubated/Non-Intubated Patients




EDAN G2 CO₂ (Sidestream)

- Superior water trap design for accurate monitoring
- iCARB™ algorithm with intelligent CO₂ pseudo wave identification technology
- Multiple sampling accessories as options for adult, child and neonate patients

Respironics CO₂ (Mainstream/Sidestream)

- Plug & play module design
- Dehumidification tube instead of water trap
- Low sampling rate of 50ml/min suitable for all types of patients

MASIMO Anesthetic Gas/O₂ (Mainstream/Sidestream)



- Unique Plug & Play module for mainstream anesthetic gas monitoring
- The no moisture sampling line design suitable for all clinical applications
- A stable low flow of 50ml/min for all types of patients



Nurse Call



VGA Output



Defib Proof & Sync



ESU Proof



Bed Rail Hook



Wall Mount



Rolling Stand

Built-in Non-volatile Memory:

A single piece of patient data contains:

1200
NIBP Measurement Review

120 h
Trend Review

200
Alarm Review

48 h
Full disclosure Wave



iM70 Patient Monitor

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Solic Medical

Ph: 416-740-8885
info@solicmedical.com



A world of potential