

Redefining Value for General Imaging Ultrasound

Acclarix™ AX4 Diagnostic Ultrasound System






Redefining innovation through value and performance

Introducing a new price performance standard for ultrasound. The remarkable Acclarix™ AX4 Compact Ultrasound System delivers a powerhouse combination of imaging performance and ease-of-use to meet the specific demands of general imaging ultrasound. The Acclarix AX4 has been designed from the ground up with a relentless focus on delivering unexpected levels of performance at a very attractive price point.


- Definitive Imaging
- Distinctive Design
- Intelligent Workflow
- Intrinsic Quality



This fluid collection on the anterolateral aspect of the knee is clearly defined with the L12-5Q transducer.



TAI for color Doppler enhances flow sensitivity and spatial resolution as demonstrated with this thyroid.



The C5-2Q transducer delivers exceptional color and spectral Doppler sensitivity visualizing the left hepatic vein.



The 3D/4D feature enables fast and easy acquisition for outstanding fetal face detail.

Stunning Clarity

Definitive imaging

Sleek design and user-friendly features are just part of the story. The Acclarix AX4 is designed to deliver amazing clarity and stunning image quality quickly—eliminating the need to tweak imaging controls.

- High fidelity architecture results in superb resolution
- Tissue Adaptive Imaging (TAI) continuously and automatically optimizes imaging allowing more focus on the patient
 - In B-mode, TAI fine tunes multiple parameters to provide the best possible image quality
 - In Doppler, TAI automatically adjusts for flow state providing improved continuity, border detection and fill-in




Redefining innovation through value and performance

Introducing a new price performance standard for ultrasound. The remarkable Acclarix™ AX4 Compact Ultrasound System delivers a powerhouse combination of imaging performance and ease-of-use to meet the specific demands of general imaging ultrasound. The Acclarix AX4 has been designed from the ground up with a relentless focus on delivering unexpected levels of performance at a very attractive price point.

- Definitive Imaging
- Intelligent Workflow
- Distinctive Design
- Intrinsic Quality

Sleek compact design fits into any environment



Tilt-and-swivel monitor for optimal viewing



Customizable control panel to personalize your user interface



Completely sealed control panel to simplify clean-up



A Clear Vision

Ultrasound reimaged for ease of use

Born of a vision to deliver meaningful design innovations that benefit the user, the Acclarix AX4 features a host of design breakthroughs that make day-to-day operation in general imaging environments easy, fast and intuitive. The result is an elegant simplicity where form and function meet at the tips of your fingers.

- Best-in-class image quality with one-key optimization
- Application-specific exam and imaging presets
- Unique gesture-control user interface
- Tilt-and-swivel HD monitor
- Completely sealed control panel aids in maintaining infection control
- Needle visualization technology improves needle identification

Transducers and Options

L12-5Q



5 - 12 MHz
Linear Array

L17-7SQ



7 - 17 MHz
High Frequency Compact
Linear Array

C5-2Q



2 - 5 MHz
Curved Linear Array

C5-2MQ



2 - 5 MHz
3D/4D Mechanical
Curved Array

E8-4Q



4 - 8 MHz
Endocavity Curved Array

MC8-4Q



4 - 8 MHz
Micro Convex Array

P7-3Q



3 - 7 MHz
Phased Array

P5-1Q



1-5 MHz
Phased Array

MC9-3TQ



3 - 9 MHz
Micro Convex Array

MT-809



Trolley System



Multi-transducer
connector option



Sturdy, versatile
travel case



Solic Medical

Ph: 416-740-8885
info@solicmedical.com



A world of potential