

# HS001

# Recovery Stretcher

## Features:

- *Locking mechanism with center-controlled locking system*
- *Steel Frame with aluminum alloy side-rails*
- *Built-in IV pole*
- *Trendelenburg: 20°*
- *Reverse Trendelenburg: 20°*
- *Backrest max angle: 75°*
- *Footrest max angle: 40°*
- *330 lbs. capacity*
- *Mattress included; only available in black colour*

## Specifications:

- **Height: 25" - 38"**
- **Length: 84.5"**
- **Width: 31.6"**

## Shipping Dimensions

- **Height: 36"**
- **Length: 87"**
- **Width: 34"**
- **Weight: 170 kg.**



\*Comes with black mattress only