

# H10

Finger Oximeter



Outpatient



Transport



**Solic Medical**

Passion to lead



**EDAN**

A world of potential

# H10 Finger Oximeter

H10 Finger Oximeter delivers a fast, reliable oximetry in a very handy, finger-tip solution by combining the monitor and the sensor into one unit. H10 is small in size, nonetheless it provides SpO<sub>2</sub>, pulse rate measurement, bar graph and plethysmogram display.



- Compact & lightweight design
- Dual color OLED display
- SpO<sub>2</sub>, PR, Bar Graph and Plethysmogram display
- One-button operation
- 4-direction display with 6 display modes
- 10-level adjustable brightness
- Low power consumption
- Automatic power off design
- Low-battery indicator
- Standard 2 AAA batteries power supply

## Features



6 Display Modes



Pocket-size Design



One-button Operation

## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**Solic Medical**

Ph: 416-740-8885  
info@solicmedical.com



A world of potential