

M3 Vital Signs Monitor



Outpatient



Transport



Ward



Portable Design



Covidien Quick Temp



M3 Vital Signs Monitor

Proven to be a durable and efficient solution for outpatients general wards, the M3 brings effective vital signs monitoring to most ambulatory and regular healthcare.

- 5.6" high resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 14 hours of continuous work
- USB flash disk support
- HL7 support via XML files or MFM-CMS

30,000
NIBP Measurement Review

72 h
Trend Review

800
Alarm Review

- Barcode scanner support
- Bi-directional communication with MFM-CMS central monitoring system
- Compatible with PatientCare Viewer PC data management software
- Wired LAN/Wi-Fi support
- Nurse call
- Suitable for adult, pediatric and neonatal patients

Seamless Connectivity



Mounting Solutions



Wall Mount



Rolling Stand



Bed Rail Clamp

Specification

Configuration	SpO ₂	NIBP	Covidien Quick Temp
SpO ₂ only	✓	—	—
NIBP only	—	✓	—
SpO ₂ +NIBP	✓	✓	—
NIBP+Covidien Quick Temp	—	✓	✓
SpO ₂ +NIBP+ Covidien Quick Temp	✓	✓	✓



Solic Medical

Ph: 416-740-8885
info@solicmedical.com



A world of potential