


# Lighten Your Vision Within a Touch

Acclarix™ AX3  
Diagnostic Ultrasound System





# The Masterpiece for Point-of-Care Ultrasound

**Redefining innovation through value and performance**

Along with the unique dual – sockets inspiration in a 4.5kg lightweight magnesium alloy body, the remarkable Acclarix AX3 Compact Ultrasound System delivers a surprising combination of features to meet the demands of point-of-care imaging. The Acclarix AX3 has been design from the relentless focus on delivering uncompromising performance at a cost-effect price to each sonographer.

- Definitive Imaging
- Intelligent Workflow
- Distinctive Design
- Intrinsic Quality



Fetal Face




Fetal Heart




MSK



Color Doppler Spleen



TVD



Anatomic M

# Stunning Clarity

## Definitive image quality

Compacted with innovative Edan TAI technology and multiple imaging processing technologies, Acclarix AX3 could perfectly display ultrasound image in different modes, assisting sonographer to make more precise diagnosis.

- High fidelity, high-channel count architecture results in superb detail resolution, particularly at depth
- Tissue Adaptive Imaging (TAI) continuously and automatically optimizes imaging allowing more focus on the patient
  - In B-mode, TAI fine tunes multiple parameters to provide the best possible image quality
  - In Doppler, TAI automatically adjusts for flow state providing improved continuity, border detection and fill-in



# The Masterpiece for Point-of-Care Ultrasound

**Redefining innovation through value and performance**

Along with the unique dual – sockets inspiration in a 4.5kg lightweight magnesium alloy body, the remarkable Acclarix AX3 Compact Ultrasound System delivers a surprising combination of features to meet the demands of point-of-care imaging. The Acclarix AX3 has been design from the relentless focus on delivering uncompromising performance at a cost-effect price to each sonographer.

- Definitive Imaging
- Intelligent Workflow
- Distinctive Design
- Intrinsic Quality

Ultimate lightweight experience with magnesium alloy body



Unique dual-sockets & dual-batteries design empowering ultra-long working time for 2 hours



Ergonomics trolley system with semi-automatic lifting



Customizable touch screen to personalize your user interface



# A Clear Vision

## Ultrasound reimaged

Born of a vision to deliver meaningful design innovations that benefit the user, the Acclarix AX3 features a host of design breakthroughs that make day-to-day operation in point-of-care environments easy, fast and intuitive. The result is an elegant simplicity where form and function meet at the tips of your fingers.



# Transducers

C5-2Q



2 - 5 MHz  
Curved Array

C5-1Q



1 - 5 MHz  
Single Crystal  
Curved Array

MC8-4Q



4 - 8 MHz  
Micro - Convex Array

MC9-3TQ



3 - 9 MHz  
Micro - convex Array

C6-2MQ



2 - 6 MHz  
3D/4D Mechanical  
Curved Array

L12-5Q



5 - 12 MHz  
Linear Array

L17-7HQ



7 - 17 MHz  
High Frequency  
Linear Array

L17-7SQ



7 - 17 MHz  
High Frequency Compact  
Linear Array

E10-3BQ



3 - 10 MHz  
Endocavity  
Curved Array ( 200° )

E10-3HQ



3 - 10 MHz  
Endocavity  
Curved Array ( 200° )

P5-1Q



1 - 5 MHz  
Phased Array

P7-3Q



3 - 7 MHz  
Pediatric Phased Array



**Solic Medical**

Ph: 416-740-8885  
info@solicmedical.com



A world of potential