

iM70 VET / iM60 VET

Veterinary Monitor



Solic Medical

Passion to lead



A world of potential

iM70 VET / iM60 VET

Veterinary Monitor



iM70 VET

12.1" Color TFT-LCD Touch Screen



iM60 VET

10.4" Color TFT-LCD Touch Screen



No-fan Design



Pitch Tone



Thermal Recorder



Shortcut Menu



Defib Proof & Sync



ESU Proof



VGA Output



Wall Mount



Rolling Stand

Standard Parameters

- 3/5-lead ECG, HR, RESP, EDAN SpO₂, EDAN NIBP, PR, 2-TEMP

Optional Parameters

- Nellcor OxiMax™ SpO₂, SunTech NIBP, 2-IBP, EDAN G2 CO₂ (Sidestream), Respironics CO₂ (Mainstream and Sidestream), Masimo Anesthesia Gas/O₂ (iM70 VET)

Multiple Display Modes



Standard Display



Large Font



Trend View



OxyCRG



Bed view

EtCO₂ Monitoring




EDAN G2 CO₂ (Sidestream)

- Superior water trap design for accurate monitoring
- iCARB™ algorithm with intelligent CO₂ pseudo wave identification technology
- Multiple sampling accessories as options

Respironics CO₂ (Mainstream/Sidestream)

- Plug & play module design
- Dehumidification tube instead of water trap
- Low sampling rate of 50ml/min

MASIMO Anesthetic Gas/O₂ (Mainstream/Sidestream)



- Unique plug & play module for mainstream anesthetic gas monitoring
- The no moisture sampling line design suitable for all clinical application
- Low sampling rate of 50ml/min

Built-in Temporary Memory:

1200

NIBP Measurement Review

120 h

Trend Review

200

Alarm Review

120 s

Frozen Waveform





iM70 VET / iM60 VET Veterinary Monitor

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Solic Medical

Ph: 416-740-8885
info@solicmedical.com



A world of potential