

# F-F05

## Transport & Recovery Stretcher

### Features:

- Hydraulic stretcher for ICU and Emergency
- Steel frame with Aluminum alloy retracting side rails, I/V pole included
- 2 gas cylinders for easy back lift control
- Trendelenburg: 18° & Reverse Trendelenburg: 18°
- Back rest max angle: 0°-90° & Footrest max angle: 0°-40°
- Weight capacity: 220Kg/485 Lbs
- 7.8" central lock silent casters
- Fifth wheel in middle for easy maneuvering
- Storage rack at the bottom with 22lbs weight capacity
- Washable blue zippered mattress



### Specifications:

- Height Adjustment: 25" x 38"
- Length: 83.5"
- Width: 30.5"

### Shipping Specifications:

- Gross Weight: 170kg.
- Carton Size: (1pc/carton): 84.5"L x 32"W x 30.5"H