

Be extra smart,
be extra refined

X10 Patient Monitor



Solic Medical
Passion to lead




EDAN

A world of potential

Build your confidence on quality care

Based on its expandability of monitoring parameters, and stunning ultra-slim design, the X10 fulfills your primary clinical requirements in various scenarios, including Emergency Rooms, General Wards, Rehabilitation Department, Cardiac Units, and In-hospital Transfer, covering all patient types, realizing from neonatal to adult monitoring.

10.1"



Standard Parameters:

3/5-lead ECG, RESP, EDAN SpO₂,
EDAN NIBP, 2-channel TEMP

Optional Parameters:


12-lead ECG, 2-channel IBP,
EDAN G2 CO₂



Proprietary algorithms & technologies

EDAN G2 CO₂ (Sidestream)

- Superior water trap design for accurate monitoring.
- iCARB® algorithm with intelligent CO₂ pseudo wave identification technology.
- Automatic zero calibration within seconds.
- Accessories for all patient types.



ECG

- Customizable filter for better performance.
- Pacemaker detection and Defibrillator/ESU protection.
- Industry leading iSEAP™ algorithm for accurate arrhythmia detection.
- SEMIP® algorithm with 208 ECG findings over age/gender diversities.

NIBP

- Dual dust filter design makes no blockage inside and provides accurate NIBP readings.
- Unique cleaning mode for routine maintenance.
- iCUFS™ algorithm with smart deflation technology.



SpO₂

- iMAT™ algorithm with outstanding motion resistance and low perfusion resistance performance.
- Reference reading of Perfusion Index (PI) from 0 to 10 according to perfusion changes.
- Simultaneous measurements of SpO₂ and NIBP of the same limb.



Nurse Call



VGA Output



Defib Proof & Sync

Built-in Non-volatile Memory:

A single piece of patient data contains:

1200

NIBP Measurement
Review

240 h

Trend Review

200

Alarm Review

48 h

Waves



X10 Patient Monitor

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Solic Medical

Ph: 416-740-8885
info@solicmedical.com



A world of potential