

Always attend to
your blood pressure

SA-10

Ambulatory Blood Pressure Monitor



Solic Medical
Passion to lead



EDAN

A world of potential

SA-10

Ambulatory Blood Pressure Monitor

The Edan ABPM SA-10 is a 24-hour ambulatory blood pressure monitor with stylish design, intelligent software, state-of-the-art blood pressure monitoring technologies, developed to avoid white coat hypertension, masked hypertension and sustained hypertension. With the help of SA-10's powerful analysis software, clinicians will manage patients' cardiovascular risks with ease and make clinical decisions with confidence.



Self-locking structure



Memory



IP22



Wireless transmission



Lightweight design



Glass panel

Feature Highlights

- Lightweight, compact and quite for better patients compliance
- Large blood pressure measurement range from 10mmHg to 290mmHg
- Proprietary algorithm with excellent motion resistance performance
- Data protection from power outage
- Auto re-measurement on failed measurement attempt

Powerful software for comprehensive analysis

With the aid of the SA-10's powerful analysis software, the clinicians can have a better insight into the variation of patient's blood pressure within 24 hours, easily identifying the periods when the blood pressure values are above limits. The informative statistical table provides a clear presentation of the measurement summary throughout the day.



Pie Chart Analysis



Circadian Analysis

	Percentage of total beds	Average	Standard deviation	Max	Min	SD limit	Nights on average	Coefficient of variation	Warning range	Alert
High (137/90)		147	21.7	175	119	56.0%		34.0%		
Medium (120/80)		125	21.7	146	100	46.0%		20.0%		
Low (103/70)	20.0%	96	5.5	98	92	5.8%		5.8%		
High (137/90)		124	21.5	150	86	18.9%		18.9%		
Medium (120/80)		27	4.5	38	20	21.1%		21.1%		
Low (103/70)		134	26.5	174	90	22.8%		19.8%		
High (137/90)		106	21.4	137	70	20.8%		22.1%		
Medium (120/80)	70.0%	57	5.8	65	50	10.0%		10.0%		
Low (103/70)		125	21.9	149	99	20.7%		22.0%		
High (137/90)		27	4.3	37	21	22.0%		22.0%		
Medium (120/80)		138	24.3	175	90	17.5%	6.5h(Normal)	17.5%	46.7	
Low (103/70)		121	21.4	146	90	18.1%	8.0h(Normal)	21.1%	20.7	0.87
High (137/90)	20.0%	57	5.5	68	50	9.7%		9.7%		
Medium (120/80)		130	21.5	150	100	19.0%		19.0%		
Low (103/70)		27	6.2	38	20	22.8%		22.8%		

Statistical Analysis



SA-10

Ambulatory Blood Pressure Monitor

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Solic Medical

Ph: 416-740-8885
info@solicmedical.com



A world of potential