

# AYR1008

# Two-Section Power Exam Table

## Features:

- *17.5" profile in low position allows multiple-level transfers*
- *80° repositioning attainable on both sides with pneumatic backrest*
- *Chair-to-Table design takes patient from consultation to examination level*
- *Heavy duty, all welded powder coated steel frame with 2" firm density foam tops*
- *Hands-free low voltage foot control*
- *Paper roll holder*
- *500 lb weight capacity*
- *Built-in casters for easy movement*
- *2-inch padded cushion in black colour*

## Specifications:

- **Height: Max: 32"; Min: 17.5"**
- **Width: 27"**
- **Length: 78"**
- **Head Section: 31.5"**
- **Seat Section: 39.5"**

