

# YXH5D

## Stretcher Chair

### Features:

- *Unique stairway mobility enabling EMS responders to move patients up or down stairs*
- *Light weight and foldable into a compact size*
- *Four 5" wheels with brakes on the rear wheels*
- *Three restraint straps with quick release buckles for patient safety during transportation*
- *2 flexible uplift rods at the front end for easy handling and 2 foldable handles at the back*
- *Weight Capacity: 182 kg*

### Specifications:

#### **Unfolded**

- **Length: 40"**
- **Width: 40"**
- **Height: 20"**

#### **Folded:**

- **Length: 40"**
- **Width: 10"**
- **Height: 20"**
- **Net Weight: 10 kg.**

